A.F.D. Animal Facility Disinfectant





Licensed Board of Veterinary Medicine CE Provider

A Multi-species Solution

Animal Facility Concentrated Disinfectant is a broad spectrum, moderate pH, hard surface disinfectant. When used as directed, this product will deliver effective biocidal action against bacteria, fungi and viruses. Biocidal performance is attained when this product is properly diluted at 1 ounce per gallon or 1:128 (See below for additional dilutions and corresponding efficacy). Animal Facility Disinfectant can be used in Veterinary Clinics, animal life science laboratories, animal research centers, animal quarantine areas, animal breeding facilities, animal holding areas, dog/cat kennels, grooming establishments, pet animal quarters, zoos, pet shops, tack shops, operating rooms, washing areas, waiting areas, examination rooms and other animal care facilities. Animal Facility Disinfectant can be used to disinfect, clean and deodorize terrarium and small animal cages, hot rock, substrate and cage furniture, plastic terrarium ornaments, driftwood, heat caves and water dishes.

Regulatory Summary		Physical Properties				
EPA Registration No.	10324-105-86550	pH of Concentrat	e 9.5 – 11.5	Fragrance	Fresh Linen	
Hospital and General Disi	nfection	Virucida	ll against (at 1 1/4	ounce per gallon)	
This product is bactericidal according to the *AOAC Use Dilution Test Method on hard inanimate surfaces modified in the presence of 5% organic serum at 1 ounce of this product per gallon of water.			This product was evaluated at 1.25 ounces per gallon use level (879 pp quat active), in the presence of 5% serum with a 10 minute contact time and found to be effective against the following viruses on hard nonporo environmental surfaces.			
Pseudomonas aeruginosa A	Canine /	Canine Adenovirus ATCC VR-800				
Salmonella enterica ATCC	Canine I	Canine Parainfluenza Virus ATC VR-666				
Staphylococcus aureus AT		Feline Calicivirus ATCC VR-782				
			ukemia virus			
Supplemental Organisms			Feline Infectious Peritonitis virus ATCC VR-2202			
			anleukopenia virus i			
Bordetella bronchiseptica A		Feline Picornavirus ATCC VR-649				
Campylobacter jejuni ATCC		Feline Rhinotracheitis virus VR-636				
Corynebacterium ammonia		Infectious Canine Hepatitis Virus ATCC VR-293				
Enterobacter aerogenes AT		Murine Parainfluenza virus type 1 ATCC VR-105				
Enterococcus faecalis ATC	C 19433	Rabies N	/ <mark>irus</mark>			
Escherichia coli ATCC 1122						
Klebsiella pneumoniae ATC		Virucidal against (at 4 ounces per gallon)				
Listeria monocytogenes AT		This product was evaluated at 4 ounces per gallon use level (2800 ppm guat active), in the presence of 5% serum with a 10 minute contact time				
Salmonella schottmuelleri A		and found to effective against the following viruses on hard nonporous				
Salmonella typhi ATCC 653		environmental surfaces.				
Shigella dysenteriae ATCC			iental sunaces.			
Staphylococcus aureus (Me	ethicillin Resistant) (MRSA).	ATCC Canine I	Parvovirus			
33591 Staphylococcus aureus Community Associated Methicillin			Minute virus of Mice (Parvovirus) ATCC VR-1346			
1 2		Porcine Parvovirus ATCC VR-742				
Resistant (CA-MRSA)	TCC 12/10					
Streptococcus salivarius A7	CC 13419	Fungici	dal against			
		This product was evaluated at 1 ounce per gallon with a 10 minute con				

Virucidal against (at 1 ounce per gallon)

Avian Influenza A (H5N1) virus Avian influenza/Turkey/Wisconsin Virus ATCC VR-798 Avian Reovirus ATCC VR-2449 Bovine Viral Diarrhea ATCC VR-534 Canine Distemper Virus ATCC VR-64 Canine Coronavirus ATCC VR-809 Equine Arteritis virus ATCC VR-796 Hepatitis B Virus Herpes Simplex Type 1 virus ATCC VR-266 Human Coronavirus ATCC VR-740 Human Immunodeficiency Virus type (HIV-1) HTLV-IIIB Infectious Bovine Rhinotracheitis virus (IBR) ATCC VR-188 Infectious Bronchitis Virus Beaudette IB42 Infectious Laryngotracheitis Virus (LT) Strain LT-IVAX Influenza A2/Japan/305 ATCC VR-100 Newcastle disease virus Norwalk/Norovirus ATCC VR-782 Porcine Respiratory & Reproductive virus (PRRSV) Porcine Rotavirus ATCC VR-893 Pseudorabies virus ATCC VR-135 Fransmissible Gastroenteritis virus Vaccinia virus Hoffmann La Roche

*The AOAC Use Dilution method is a method of testing the efficacy of disinfectants, originally developed in 1955. Throughout its numerous revisions, it has become the standard for evaluating liquid and dilutable liquid disinfectants for hard surfaces. In particular, this method is specified by the U.S. EPA as the required method for disinfectant claim substantiation.

This product was evaluated at 1 ounce per gallon with a 10 minute contact time and found to be effective against the following fungi on hard nonporous environmental surfaces.

Trichophyton mentagrophytes ATCC #9533

Mold and Mildew Control

Use this product at 1 ounce per gallon to control the growth of mold and mildew and their odors on hard, nonporous surfaces. Thoroughly wet all treated surfaces completely. Let air-dry. Repeat application weekly or when growth or odor reappears.

Aspergillus niger ATCC #6275

Non-Food Contact Surface Sanitizer

Add 1 ounce of this product to 1 gallon of water to sanitize hard nonporous non-food contact surfaces. Treated surfaces must remain wet for 60 seconds. Then wipe with sponge, mop or cloth and allow to air dry. At this dilution, food contact surfaces must be rinsed

Klebsiella pneumoniae ATCC 4352 Staphylococcus aureus ATTC #6538

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Apply this product with a cloth, mop or mechanical spray device. When applied with a mechanical spray device, surfaces must be sprayed until thoroughly wetted. Treated surfaces must remain wet for 10 minutes. For sprayer applications, use a course pump or trigger sprayer. Spray 6-8 inches from surface. Rub with brush, sponge or cloth. Prepare a fresh solution at least daily or when use solution becomes visibly dirty. For heavily soiled areas a preliminary cleaning is required. Please refer to product label and SDS for additional product use and safety information.